# SOPIRA® Carpule® Syringe Aspirating

<b>DE</b>	Gebrauchsanweisung.....	2
<b>GB</b>	Instructions for use .....	6
<b>FR</b>	Mode d'emploi .....	10
<b>ES</b>	Instrucciones de uso.....	14
<b>IT</b>	Istruzioni per l'uso .....	18
<b>PT</b>	Instruções de uso .....	22
<b>NL</b>	Gebruiksaanwijzing.....	26
<b>SE</b>	Bruksanvisning.....	30
<b>NO</b>	Bruksanvisning.....	34
<b>GR</b>	Οδηγίες Χρήσης .....	38
<b>CZ</b>	Návod k použití.....	42
<b>PL</b>	Instrukcja obsługi .....	46
<b>HR</b>	Upute za uporabu .....	50
<b>BG</b>	Инструкции за употреба.....	54
<b>EE</b>	Kasutusjuhend .....	58
<b>SI</b>	Navodila za uporabo .....	62

**KULZER**

mitsui chemicals group

Manufacturer:  
Henke-Sass, Wolf GmbH  
Keltenstrasse 1  
78532 Tuttlingen (Germany)

Distributor:  
Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau (Germany)

Made in Germany

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Instructions for use **GB**

## Aspirating dental cartridge syringe

In combination with SOPIRA Carpule disposable needles the SOPIRA Carpule syringe provides a universal, easy-to-use injection system.

- All standardised cartridges (1.8 ml) with rubber plunger (solid as well as recessed) can be used with the SOPIRA Carpule Syringe.
- The rotatable cartridge holder permits perfect view of the cartridge in every situation, especially when controlling the aspiration process and the volume of the injected solution.
- The chromium-plated SOPIRA Carpule Syringe is easily to dismantle, allowing easy cleaning and sterilization.

## Contraindication


Use is contraindicated where there is a known or presumed allergy to components in the product.

## Side effects


The product or one of its components may cause hypersensitive reactions in individual cases. You can request information on the ingredients if you suspect this is the case.

**Warning: For use as intended by dental professionals only.**


## Instructions for use:


Withdraw the piston rod completely.

Load the cartridge.


Screw the spiral hooks clockwise into the rubber plunger applying light pressure.

**Fix the needle:**

Ensure that the cartridge end of the needle pierces the centre of the rubber membrane. It is advisable to insert the needle just prior to injecting. Expel a few drops of the solution. Before injecting, check that the aspiration help is firmly and correctly positioned in the barrel used. The syringe is now ready for use.

**After injection:**

Detach the needle and remove the spiral hooks from the rubber plunger. Fully retract the piston rod and remove the cartridge.

## Preparation:

The syringe must be prepared prior to initial use and each subsequent use (syringes are supplied clean but not sterile):

<b>Preparing at point of use:</b>	The syringe must be checked for any damage such as coating each time before use. Syringes must no longer be used if they are damaged. Remove any coarse dirt from instruments directly after use. Do not use any fixing agents or hot water (>40 °C) as this causes residues to deposit and can have an impact on whether cleaning is successful or not.
<b>Transport:</b>	Safe storage and suitable transport in a closed container are recommended to prevent damage to the medical product and avoid any contamination of the environment.
<b>Manual pre-cleaning:</b>	Use a soft, multi-purpose brush to clean the medical products under a cold running tap until all visible residues and contaminants are removed. Rinse any holes, threads and cavities for at least 10 sec. with a water pistol at a pressure of 3.8 bar.
<b>Cleaning:</b>	Place medical products on a screen tray on a loading rack and start the cleaning process <ol style="list-style-type: none"><li>1. Pre-wash with cold water for 4 minutes</li><li>2. Drain</li><li>3. Pre-wash with 0.5 % Neodisher Mediclean, Dr. Weigert (Hamburg) for 10 min.</li><li>4. Drain</li><li>5. Neutralize with 0.1 % Neodisher Z, Dr. Weigert (Hamburg) (&gt; 40 °C) for 6 min.</li><li>6. Drain</li><li>7. Rinse with demineralized water (&gt; 40 °C) for 3 min.</li><li>8. Drain</li></ol>
<b>Disinfection:</b>	Carry out machine thermal disinfection while taking into account national requirements regarding the AO value (see ISO 15883).
<b>Drying:</b>	The washer-disinfector drying cycle is used to dry the exterior of medical products. A lint-free cloth can be used to dry the product further by hand. Sterile compressed air must be used to dry cavities in medical products
<b>Functional testing and maintenance:</b>	A visual inspection must be conducted to ensure that the product is clean. It must then undergo maintenance and a functional test as per the operating instructions. If necessary, reconditioning must be repeated until the medical product is visibly clean.
<b>Packaging:</b>	Standard-compliant packaging for medical products used for sterilization in line with ISO 11607 and EN 868
<b>Sterilization (autoclaving):</b>	Products must be sterilized using a fractionated pre-vacuum procedure (as per ISO 13060/ISO 17665) while taking into account requirements applicable in the country of use. <ol style="list-style-type: none"><li>1. Fractionated pre-vacuum (triple)</li><li>2. Sterilization temperature of 134 °C</li><li>3. Shortest exposure time: 3 min. (full cycle)</li><li>4. Drying time: at least 10 min.</li></ol>
<b>Storage:</b>	Sterilized medical products must be stored in a dry, clean, dust-free environment at moderate temperatures between +5 °C and +40 °C.
<b>Information on the validation of preparation:</b>	The following machines and materials were used for the validation: Cleaning agent: Neodisher Mediclean (alkaline); Dr. Weigert; Hamburg Neutralizing agent: Neodisher Z; Dr. Weigert; Hamburg Washer-disinfector: Miele G 7836 CD

<b>Additional instructions:</b>	According to the German Act on Medical Devices, users are responsible for the validation of their preparation processes. This also applies if the chemicals and machines specified above are not available. The user must ensure that the reprocessing method, including resources, materials and personnel, are suitable for achieving the required results. The state of the art and national laws require the observance of validated processes.
---------------------------------	---

Do not undo or remove any screws (tilt mechanism/finger ring) before preparation.

According to the EU Medical Device Regulation, users/patients are obliged to report any serious incidents with a medical product to the manufacturer and the relevant authorities in the country where the incident occurred. Observe the safety data sheet and national regulations when disposing of the device. Always indicate the batch number and item code when reporting incidents involving the product.

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