dima[®] Print Denture Teeth

- DE** Gebrauchsanweisung
- GB** Instructions for use
- FR** Mode d'emploi
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- PT** Instruções de uso
- NL** Gebruiksaanwijzing

Explanation of symbols on labeling

- LOT** Batch code
- Manufacturer
- Keep away from sunlight
- Consult instructions for use
- Use-by date
- Storage temperature

Manufacturer:
DENTCA, INC.
357 VAN NESS WAY, STE 250
TORRANCE, CA 90501
USA

Distributor:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)
+49 800 4372522

Distributed in USA / Canada exclusively by:
Kulzer, LLC
4315 South Lafayette Blvd.
South Bend, IN 46614-2517
1-800-431-1785

Caution: Federal law restricts this device to sale by or on the order of a dental professional.

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EU Representative:
MT Promed Consulting GmbH
Altenhofer Str. 80
66386 St. Ingbert
Germany

Medizinprodukt
Medical device
Dispositif médical
Producto sanitario
Dispositivo medico
Equipamento para saúde
Medisch hulpmiddel

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GB Instructions for use dima[®] Print Denture Teeth

Indications for Use
dima Print Denture Teeth is a light-curable polymerizable resin to fabricate, by additive manufacturing, preformed denture teeth to be used in a denture. The fabricated tooth is an alternative to preformed plastic tooth for denture. The fabrication of denture tooth with dima Print Denture Teeth requires digital denture tooth files instead of physical moulds, a stereolithographic additive printer, and curing light equipment.

Requirements

Digital denture Teeth; STL format file	Operation Software	Provider
Additional printer and its operation software; <td>cara Print CAM</td> <td>Kulzer</td>	cara Print CAM	Kulzer

Curing light equipment	Provider
HiLite power 3D	Kulzer

Specific Manufacturing Considerations

- Digital denture tooth file
 - File format: STL file
 - Digital design: DENTCA BYTE tooth – Universal, Oval, and Square
 - File size: Upload-able in the 3D printer operation software
- Additional printer
 - Hardware
 - Wavelength : 385 nm or 405 nm
 - Light source
 - Stereolithographic (SLA) method; laser with 25 mW < X < 250 mW
 - Digital Light Processing (DLP) method; high power LED or lasers
 - Build Volume: > 103 x 58 x 130 mm (Least fit one arch)
 - Build Path: line drawing path or surface layer drawing path
 - Features of Operation Software
 - STL file import
 - Automatic rotation and placement
 - Layer slicer for path inspection
 - Auto and manual generation of supports
 - Printing Parameter

Printer Model	Layer Thickness (micron)	Recommended orientation angle (degree)	Support point size (mm)	Support density (mm)
cara Print	50, 100	20-40	1-2	min. 1

- Environmental Conditions
 - Temperature: 18–30°C (64–86°F)
 - Relative Humidity: 30–90%
- Cleaning Kit
Rinse bath and tubs, flush cutter, paper towel, squeeze bottle for isopropyl alcohol, Scraper
- Recommended Curing light equipment (Post curing units)
 - Flood Type Curing Equipment

Manufacturer / Model	Curing Chamber	Supply voltage (V / Hz)	Lamp power (W)	Light intensity (mW/cm²)	Lamp wavelength (nm)	Curing time (min)
Kulzer / HiLite power 3D	Required	100,115,230 / 50–60	200	N/A	390–540	20

- Accessories
 - USP Grade glycerin
 - Transparent glass container and 2 glass plates
 - Heat-protective gloves and silicone coated stainless steel tong
 - Thermocouple
- Notification
The device specifications have been validated using the software, printers, and process parameters specified in this document. Any other printers, operation software and post-printing processes will be outside of the device specifications and the FDA clearance. Users shall follow this document to use the device.

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Warnings:
1. dima Print Denture Teeth contains polymerizable monomers which may cause skin irritation (allergic contact dermatitis) or other allergic reactions in susceptible persons. If contact with skin, wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance.
2. Avoid inhalation or ingestion. High vapor concentration can cause headache, irritation of eyes or respiratory system. Direct contact with eyes may cause possible corneal damage. Long-term excessive exposure to the fumes may cause more serious health effects. Meet or exceed OSHA standards. Eye Contact: Immediately flush eyes with plenty of clean water for at least 20 minutes, and consult a physician. Wash the contact area thoroughly with soap and water. Inhalation: In case of exposure to a high concentration of vapor or mist, remove person to fresh air. Give oxygen or artificial respiration as required. Ingestion: Contact your regional poison control center immediately.
BURN HAZARD: GLYCEROL BATH CAN REACH TEMPERATURES OF 90°C (~200°F) AND LEAD TO SEVERE BURNS. Only trained users should perform the glycerol curing step with caution and appropriate PPE. We also recommend placing a warning label on the window of the cure unit to alert all lab users to the potential hazard.

Precautions:

- When washing the printed denture Teeth with solvent or grinding the denture Teeth, it should be in a properly ventilated environment with proper protective masks and gloves.
- Store dima Print Denture Teeth at or below 15–25°C (60–77°F) and avoid direct sunlight. Keep container closed when it is not in use. Product shall not be used after expiration date.
- Expired or unused dima Print Denture Teeth should be completely cured or polymerized prior disposal.

Adverse Reactions:

- Direct contact with the uncured resin may induce skin sensitization in susceptible individuals.
- Proper ventilation and personal protective equipment should be used when grinding printed denture Teeth as the particulate generated during grinding may cause respiratory, skin and eye irritation.

Procedure to Fabricate the Denture Teeth

- Printing Preparation
 - Select the denture Teeth shade based on prescription. (Recommended to use the different resin tank or tray for the different shade.)
 - Open the 3D printer cover and fill the resin tank or tray of the printer with dima Print Denture Teeth up to the required filling line by the manufacturer. (When filling the resin into the resin tank or tray, gloves and mask should be used.)
 - Close the printer cover.
 - Printing
 - Load the denture Teeth model file to be printed in printer operation software which printer manufacturer recommended.
 - Use the software tool to rotate the model in order to locate the model in proper position to the build plate form.
 - Tilt the lingual side and root of Teeth to around 20 to 40 degree against the build plate form.
 - Generate support sticks on the denture Teeth using the recommended setting by printer provider. The support is not enough, add supports on the Teeth. (Avoid the support structures on the valley between Teeth.)
 - Start printing.
 - Cleaning
 - Detach the printed Teeth from the build platform.
 - Use a small flush cutter to remove the support sticks from the denture Teeth.
 - Wash the denture Teeth with isopropyl alcohol.
 - Use air blowing to dry the denture Teeth or dry it at room temperature under ventilation system or open area.
 - Teeth Post Curing
 - Smooth the support marks using a bur after washing the denture Teeth with water and drying.
 - For printed denture fabrication, skip this post-curing step and move to Denture Fabrication using a printed Teeth and "base" step. For conventional fabrication to use as preformed Teeth, the printed Teeth should be cured by soaking into glycerin container for the required curing time under recommended post-curing unit.
 - Use the printed preformed Teeth to fabricate a denture by conventional denture process.
- Denture Fabrication using a printed denture Teeth and base
 - Place the printed Teeth to the printed denture base.
 - Prepare, before post-cure, the printed denture base with socket shapes to receive printed Teeth (Tooth Sockets).
 - bonding the printed Teeth into the corresponding Tooth Sockets on the printed denture base and check Teeth fitting.
 - Apply the small amount of light curable adhesive into the Tooth Sockets and bond Teeth by exposing into UV light until the Teeth set in position. Recommended to use dima Print Denture Base for bonding.
 - If necessary, apply small amount of dima Print Denture Base using an applicator to smoothen the edges of the denture base and cure it.
- Post Curing
 - Cure the final denture by sinking into the glycerin container (glycerin temperature should be greater than 60°C and it is recommended to replace every 80 hrs or every three months whenever comes first or until the required curing time under recommended post-curing unit. For half of the post-curing time issue side up and for another half of the time issue side down.
 - Take out the printed denture from curing oven using coated tong. (Be careful hot glycerin!)
 - Rinse the cured denture with a water.
- Finishing
 - Smoothen the support spots on the denture base using a bur.
 - Polish the final denture with wet polishing sand by conventional method.

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