PALA®

PalaVeneer® Dentine

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GB Instructions for use

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SE Bruksanvisning

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NO Bruksanvisning

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KULZER  
MITSUI CHEMICALS GROUP

Distributed by:  
Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau (Germany)  
+49 800 4372552

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Stand: 2021-07

99001600000

## GB Instructions for use PalaVeneer® Dentine

Product description  
Tooth colored acrylic for prosthetics, self curing.  
PalaVeneer Dentine is a self-curing, tooth colored PMMA prosthetic acrylic, colour-matched to Kulzer PalaVeneer veneer shells, for all indications listed below. It consists of a powder component in various shades (A1 – D4, 2 bleach shades, 1 incisal material) and a liquid component.  
Prosthetic acrylic for casting permanent and removable dental restorations.

Composition  
Main powder component:  
methyl methacrylate copolymer.  
Main liquid components:  
methyl methacrylate, dimethacrylate.  
PalaVeneer Dentine is cadmium free.

Indications  
• Fixing and backing PalaVeneer veneer shells on metal frameworks  
• Fixing and backing milled acrylic teeth on metal frameworks  
• To restore/bond facets  
• To produce long term temporaries (facet technique)  
• To aesthetically finish pre-fabricated teeth

Contraindications  
Patient contact with unpolymersed PalaVeneer Dentine is contraindicated. If the patient is allergic to one or more ingredients of PalaVeneer Dentine, the product must not be used.

Advice  
Only for processing in dental laboratories in accordance with intended use. The product contains polymerisable monomers (e.g. methacrylates, acrylates), which may sensitise the skin of sensitive people. If any irritation appears or there are known allergies to resins based on methacrylates or acrylates, the product should not be used.

Fixing and backing PalaVeneer veneer shells/milled acrylic teeth/facets on metal frameworks

Setup  
Adapt the PalaVeneer veneer shells to the framework structure. Set-up the veneer shells with soft, dentine-shade wax, such that position corrections are still possible. On completion of the setup, fill up and model cervical, approximal, palatal/lingual cavities with hard wax.

Matrix production  
Create the matrix from kneadable impression silicone with increased hardness of Shore A > 60. Remove the silicone matrix after curing.

Conditioning the metal framework  
We recommend mechanical retention on the metal framework. Clean and blast the metal

framework (110 – 125 µm Al<sub>2</sub>O<sub>3</sub>, 3 bar), blow off dust with oil-free compressed air, do not steam off. Condition metal framework with Signum® metal bond I and II (see Signum® metal bond I and II instructions for use). Coat metal framework with 2– 3 thin layers of light-curing Signum® opaque F and polymerise each layer individually. For an optimum bond, apply the opaque without undue delay and without contaminating the surface by touching it. Then apply the PalaVeneer Dentine.

Conditioning  
Remove veneer shells, steam off wax residues from veneer shells, blast the contact surfaces of the veneer shells to the PalaVeneer Dentine (110 – 125 µm Al<sub>2</sub>O<sub>3</sub>, 2 bar) or roughen with a coarse diamond. Blow off dust with oil-free compressed air, do not steam off. Fix the veneer shells in the silicone matrix.  
To improve the bond between PalaVeneer Dentine and the veneer shells, all contact zones of the veneer shells to the PalaVeneer Dentine are wetted with the Palabond® adhesion primer. Apply Palabond® twice using a brush (without metal ferrule) and allow to set for 30 seconds each time. After the second coat, the adhesion primer remains active for approx 10 minutes.

Preparing the models  
After removal of the wax with hot water (without chemical additives), the stone model is isolated twice with thin layers of Aislar® (see Aislar® instructions for use). This is the only way to ensure perfect surface isolation between the model and PalaVeneer Dentine.

Dosing and processing  
PalaVeneer Dentine is processed with a mixing ratio of 1.8 g powder : 1 g liquid. Introduce the liquid to the mixing beaker, add the appropriate quantity of powder rapidly and stir for approx. 30 seconds to make a homogeneous paste. Prevent bubbles from forming.  
PalaVeneer Dentine can be processed after mixing at a room temperature of 23°C (73°F) for approx. 3 min. The processing time depends on ambient temperature.

Pre-pressing technology: Apply PalaVeneer Dentine in the fixed veneer shells with a slight excess, then pre-press and fix the silicone matrix. Remove excess and model the palatal/lingual surfaces.

Polymerisation  
Polymerisation time in the Palamat elite polymerisation unit: 20 min. water temperature: 55°C (131°F), pressure: 2 bar.

Trimming and polishing  
PalaVeneer Dentine is trimmed e.g. by means of a carbide bur, lens-shaped silicone burber polisher, separating disc, goat-hair bristles with plastic polishing paste and a cotton buff with high-gloss polishing paste and a speed of 5,000 – 10,000 rpm.

Aesthetic finishing of prosthetic teeth  
Roughen pre-fabricated teeth to be finished using a medium-fine diamond or by blasting (110 – 125 µm Al<sub>2</sub>O<sub>3</sub> and 2 bar pressure) and condition twice with Palabond liquid. Mix PalaVeneer Dentine then build up and polymerise the material at the relevant points.

Storage advice  
Do not use the material after the expiration date. After opening, do not store above 25°C (77°F). Avoid direct sunlight. Close containers immediately after use. Keep out of reach of children. Before using the product for the first time, ensure that the primary packaging is unbroken. In case of damage, do not use.

Hazard/Safety information liquid (methyl methacrylate, tetramethylene dimethacrylate)

Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat/sparks/open flames/hot surfaces. – No smoking. Wear protective gloves/protective clothing/eye protection/face protection. Ground/bond container and receiving equipment. Take precautionary measures against static discharge.

When reporting back to the product, please always state the batch number and article number. According to the EU Medical Device Regulation, users / patients are obliged to report serious events with a medical device to the manufacturer and the competent authority in the country in which they occurred.

Recommendations  
Disposal according to official regulations. Do not dispose of contents and packaging that has not been completely emptied with household waste or allow them to enter the sewage system. European waste catalog: 180106 Chemicals consisting of or containing dangerous substances.

Safety data sheets and more information are available at our website www.kulzer.com

© = registered trademark of Kulzer GmbH  
Dated: 2021-07

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© = geleverd onder handelsmerk van Kulzer GmbH
© = indregistreret varemærke tilhørende Kulzer GmbH
© = rekisteröity tavaramerkki Kulzer GmbH