



Venus®
Venus® Flow
Simply beautiful. Simply Venus.

Venus helps you give your patients long-lasting, aesthetically pleasing restorations. With this universal composite you can create durable, gleaming restorations using the simplest techniques. Venus Flow adds a perfect flowable variant to the outstanding properties of the Venus range. With this well-attuned system you are ready to face the greatest challenges posed by anterior and posterior restorations. Its comprehensive range of shades and unique layering concept will ensure that all your restorations are masterpieces. Why should you settle for less?

Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

Venus® – Simply beautiful. Simply Venus.


Restorations with no ifs or buts

The refractive indices of the fillers and matrices in our Venus products have been perfectly aligned. Our 'Color Adaptive Matrix' is responsible for this. It gives Venus the unique ability to blend in with the surrounding tooth structure, creating an outstandingly natural look. The result of this complex mechanism is virtually undetectable restorative margins. The ultrafine filling particles allow you to achieve the perfect polish comparable to natural enamel.

Universal application

Venus is a medium-viscosity filling material with exceptional handling properties. The Venus composite is stable whilst at the same time being smooth to work with. This facilitates contouring and increases efficiency. The excellent polishing and layering properties of Venus make it ideal for anterior restorations. Its extremely good strength and abrasion resistance also make it the material of choice for posterior fillings.

The classic universal composite: evidence-based and clinically tested



Source: Dunn J: 2 year clinical evaluation of venus composite in anterior teeth. J Dent Res 86 (spec Iss A), 1575, 2007.

	Venus			Venus Flow
	1 x 4g Syringe	20 x 0.25g PLT	10 x 0.25g PLT	1 x 1.8g Syringe
A1	X	X		X
A2	X	X		X
A3	X	X		X
A3.5	X	X		X
A4	X		X	X
B1	X	X		
B2	X	X		X
B3	X		X	X
C2	X	X		
C3	X		X	
C4	X		X	
D2	X		X	
D3	X		X	
HKA2.5	X	X		X
HKA5	X		X	
OA2	X	X		X
OA3	X		X	
OA3.5	X		X	
OB2	X		X	
T1	X		X	
T2	X		X	X
SB1	X		X	X
SBO	X		X	
Baseliner				X

Venus® Flow – Low viscosity, high performance.

Venus Flow is the perfect addition to Venus composite in clinical situations in which precise control can only be attained by using a flowable composite. With this well-attuned system you are ready to face the daily challenges posed by anterior and posterior restorations.

Venus Flow is easy to use and, thanks to its radiopacity, ensures a reliable diagnosis. Its low viscosity and good thixotropic properties make Venus Flow the material of choice in a diverse range of indications, such as sealing fissures, lining cavities and creating smaller fillings.

