

**Delara®**

GB Mould chart  
DE Formenkarte  
FR Cartes des Formes  
IT Tabella forme  
ES Carta de Formas  
PT Carta Molde  
NL Vormenkaart  
SE Formkort  
DK Formkort  
FI Muotto korti  
NO Formkort  
EE Vormikaart  
LV Formas karte  
LT Formos kortelė  
ARA مخطط القالب



**Delara®**

A modern all-round tooth line for all indications with lifelike aesthetics.

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP

## Everybody's Delara



Delara meets every need offering a lively surface structure and unique layering, along with modern and bright shading for a fresh and vivid appearance.

- **Easy handling for predictable results:** freedom during set-up without excessive grinding or changing the basic shape and layering effects. The polishing buffer, basal grinding reserve and optimized incisal edge simplify the set-up process.
- **Flexibility:** bigger interdental space, full range of shades, process-optimized and anatomical forms – for all indications and set-up options.
- **Aesthetics:** blueish effect for a fresh and lively appearance, along with modern shading, a vivid surface structure and dynamic light refraction.

That's why we call it  
"Everybody's Delara".

*“Delara anterior teeth look so natural and beautiful in the mouth thanks to the unique surface structure. In addition, they are easy to set up, especially the posterior teeth. With Delara you have the freedom of set-up with almost no grinding necessary.”*



**Ilkka Garaisi**

As dentist (clinical dental technician) he has over 25 years of experience and treated over 7.000 patients with full dentures. He is president of EHT-federation in Finland, teacher of denturists, and representative of Finland in the IFD.



# Anteriors

Frontzähne / Antérieures / Anteriori / Anteriores / Anteriores / Tanden / Framtänder / For tænder / Etuhampaat / Fortennene / Eesmised hambad / Priekšējie zobi / Priekiniai dantys / الأمامية

<b>T 42</b>   	42,4 9,6 8,0	<b>S 44</b>   	44,5 10,4 8,3
<b>O 43</b>   	44,2 10,5 8,3	<b>T 44</b>   	45,1 10,5 8,2
<b>S 43</b>   	43,9 9,6 8,2	<b>R 45</b>   	45,0 10,6 8,2
<b>T 45</b>   	45,7 11,3 8,4	<b>O 47</b>   	48,1 11,4 9,1
<b>T 46</b>   	46,5 10,7 8,7	<b>R 47</b>   	47,5 11,5 8,8
<b>R 46</b>   	46,5 10,1 8,8	<b>S 48</b>   	49,0 11,3 9,1
<b>R 48</b>   	48,8 13,0 9,3	<b>S 50</b>   	51,0 11,8 9,6
<b>T 49</b>   	49,6 12,2 9,1	<b>O 51</b>   	52,2 12,4 9,8
<b>R 49</b>   	48,7 10,9 9,3	<b>R 52</b>   	52,1 13,6 9,8



## Posteriors

Seitenzähne / Postérieures / Posteriori / Posteriores / Posteriores / Kiezen / Bakre tänder / Bageste tænder / Takimmaiset hampaat / Bakre tenner / Tagumised hambad / Aizmugurējie zobi / Užpakaliniai dantys / الخلفية





Upper anteriors

Lower anteriors

Posteriors



T 42
T 44
T 45
T 46
T 49

L 33
L 33 + L 35
L 33 + L 35 + L 36
L 35 + L 36 + L 38
L 38 + L 39

30 + 32
30 + 32
32 + 34
32 + 34
36



O 43
O 47
O 51

L 33 + L 35 + L 36
L 38 + L 39
L 38 + L 39

30 + 32
34
36



R 45
R 46
R 47
R 48
R 49
R 52

L 33 + L 35
L 33 + L 35 + L 36
L 38 + L 39
L 36 + L 38
L 38 + L 39
L 38 + L 39


32
32 + 34
34
34
34 + 36
36




S 43
S 44
S 48
S 50

L 33 + L 35 + L 36
L 33 + L 35 + L 36 + L 38
L 38 + L 39
L 38 + L 39

30 + 32
30 + 32 + 34
34 + 36
36




**T 42**  
42,4  
9,6  
8,0



**30 U**  
29,1  
9,1



**L 33**  
33,6  
9,4  
5,2



**30 L**  
31,4  
8,1

This mould chart complies with EN ISO 22112:2017. / Diese Formenkarte erfüllt EN ISO 22112:2017. / Cette carte forme remplit EN ISO 22112:2017. / Questa tabella forme soddisfa EN ISO 22112:2017. / Esta carta de formas es conforme con la EN ISO 22112:2017. / Esta carta molde está em conformidade com EN ISO 22112:2017. / Deze vormenkaart voldoet aan EN ISO 22112:2017. / Detta formkort uppfyller EN ISO 22112:2017. / Dette formkort oppfyller EN ISO 22112:2017. / Tämä muotokortti täyttää EN ISO 22112:2017. / Dette formkortet oppfyller EN ISO 22112:2017. / See vormikaart vastab standardile EN ISO 22112:2017. / Ší formos karte atbilst EN ISO 22112:2017. / Ši formos korta išpildo EN ISO 22112:2017. / يتوافق مخطط القوالب هذا مع EN ISO 22112:2017.

**Manufacturer**

Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau, Germany

**Contact in Germany**

Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau, Germany  
info.lab@kulzer-dental.com